

Miniball electronics at CERN - May 2005

Contents

1	Foreword	3
2	Notation	3
3	The delayed CD quadrant signal	3
4	Generation of the gamma gate	4
5	Generation of triggers	4
6	Generation of the DAQ dead signals	5
7	Generation of the DAQ trigger	5
8	Generation of the ADC gates and the TDC stops	6
9	Generation of the EBIS window	6
10	Generation of the off window	7
11	Generation of the on/off window	8
12	The generation of the GFLT	8
13	The DGF BUSY/SYNCH loop	9
14	Note about special signals	9
15	Generation of the T1 timestamp	9
16	Generation of the proton supercycle timestamp	10
17	Generation of the control signals for the scalers	11
18	Generation of the control signals for the pattern unit	11
19	Generation of the forced readout	11
20	The pattern unit bits	12
21	Laser on/off	12
22	1 MHz and EBIS, 1 MHz and GFLT	12
23	Ionisation chamber	13
24	Scalers	14
25	Positions of modules in crates and racks	14
25.1	Rack 1	14
25.1.1	CAMAC crate R1.C2	14
25.1.2	NIM crate R1.C3	14
25.2	Rack 2	15
25.2.1	KM-6 crate R2.C1	15
25.2.2	KM-6 crate R2.C2	15
25.2.3	NIM crate R2.C3	15
25.3	Rack 3	16
25.3.1	VME crate R3.C1	16
25.3.2	NIM crate R3.C2	16
25.3.3	NIM crate R3.C3	17

¹Nigel Warr May 2005

25.3.4	NIM crate R3.C4	17
25.4	Rack 4	17
25.4.1	CAMAC crate R4.C1	17
25.4.2	CAMAC crate R4.C2	18
25.4.3	R4.C3	19
25.4.4	NIM crate R4.C4	19
25.5	Rack 5	19

1 Foreword

This document describes the state of the electronics used at CERN for the May 2005 setup. It is based on the one for July 2004, though in a few places changes made for subsequent experiments have been indicated, notably the addition of laser on/off bits in the pattern unit and the Munich ionisation chamber. That means it is different to all previous experiments and is likely to be different to future ones. However, there should be certain similarities.

The electronics for this setup has grown gradually rather than being thought out as a single setup, so modules were not necessarily placed in the optimal position in the crates. In places, signals were duplicated due to different people building different parts of the setup and not realizing that the signal they were constructing was already available elsewhere.

In November 2004, a part of the electronics was dismantled and sent to GSI, Darmstadt, and the racks were all moved in order to allow for the building work to extend the ISOLDE hall. However, care was taken then to label things, so that they could be put back together. In May 2005, we put the things back together and it is believed that the state of November 2004 was achieved. In doing so, we noticed some changes with respect to the July 2004 setup. These all seem to be deliberate improvements.

However, the setup works, so there is no reason to change it for this campaign, but it could be improved on in the future.

Caveat: this document was produced during the setup with all the interruptions that that entails, so it probably contains mistakes. Several mistakes have already been found and corrected. There are likely to be others! Also people were still modifying the electronics while I was noting down the configuration. Use it with caution!

The main difference between the 2003 and 2004 setups is that in 2003 all of the CD was treated as one unit, whereas in 2004, each quadrant was dealt with separately. Also, in 2004 we took the ADC busy directly from the ADC module, whereas in 2003 we generated our own signal which was longer than the one generated by the module. Back then, we thought (erroneously) that the ADC was giving a busy signal which was too short and was the reason for mismatches in the number of ADC gates and the number of timestamper DGF events. In fact, the problem turned out to be due to the DGF sometimes missing the last event, so there was no reason to artificially generate the ADC busy.

2 Notation

I have used the notation “R2.C4.S13” to mean rack two (starting at 1 nearest the wall) crate four (starting with 1 at the top), slot 13. Similarly “R2.C4” means all of crate four of rack two. In the diagrams, the location of the modules is indicated in magenta.

3 The delayed CD quadrant signal




Figure 1: The generation of the delayed CD quadrant signal. This is done for each of the four quadrants of the CD.

From the CD we need to generate logic signals to indicate when an event has occurred in a given quadrant. Note that in previous setups we treated the whole CD as one, but now we treat it as four independent quadrants.

For each quadrant, we take the logic output from the RAL 109 shaper amplifiers via a flat cable to four CAMAC MALU units. Note that although these are CAMAC modules, we do not need to program them, so they are in a crate without a CAMAC crate controller. They just take their power from the CAMAC crate.

Each MALU gives a differential ECL ORO output (this is just an OR of all the inputs) and these signals for all four MALUs are combined on a single flat cable and sent to an ECL → NIM converter. From there,

the signal for each quadrant goes into one quarter of a 4x4 fan-in/fan-out module and then is used to start a dual timer.

We use three copies of the signal from this fan-out.

- One goes to a scaler. We use scaler 65 for quadrant 1, scaler 66 for quadrant 2 etc.
- One copy is delayed by 800 ns to generate the delayed CD quadrant signal.
- The third copy is converted to ECL and used to start the TDC corresponding to that quadrant.

We want to delay this signal by 800 ns, but the CAEN dual timers have the property that if a second signal arrives during that time, the end marker comes 800 ns after the **second** signal, not the first. So we use half of a dual timer to generate a signal which is 800 ns wide and then use that to start the second half of the same module also set to 800 ns (so both start more or less simultaneously, but the delaying of the end marker is suppressed). We then take the end marker from the second half of the module. There are probably better ways to do this. In all probability this behaviour of the dual gate is a feature but it doesn't seem to be possible to turn it off and the manual is utterly useless. Perhaps a different kind of gate generator would be better, or a discriminator followed by a gate generator. In future years we should really rethink this part.

For all four quadrants, we use four MALUs, one ECL → NIM converter, one fan-in/fan-out module and four dual timers.

Note that the RAL109s and the NIM crate containing the NIM → ECL converter are in the CD electronics rack, and the CAMAC crate with the MALUs are in the leftmost rack.

Note that a special cable is used from the NIM → ECL unit to the TDCs which takes the four signals (one from each quadrant) and split so that each one goes to start a different CAEN V775 TDC.

4 Generation of the gamma gate




Figure 2: The generation of the gamma gate. This is the logical OR of all the signals from the Ge cores.

We need to generate a signal indicating that one of the Ge detectors has an event. To do this we take the Mult Out (35 mV per hit) from the DGF having the first core signal and put it into the Mult In of the next DGF with a core signal (i.e. the third DGF) and from its Mult Out to the DGF with the third core. We then send this signal (i.e. one for each of the eight clusters) to a specially modified TFA (which contains eight Ge preamplifiers set to saturate with a 35 millivolt signal and generate a -0.8 Volt NIM output). In other words, the special module simply converts the 35 mV per hit signal into NIM logic. This signal is passed into a fan-in/fan-out.

This signal is used to generate the particle-gamma coincidence. Note that its width is determined by the DGF paramter *FTPWIDTH*.

5 Generation of triggers

The delayed CD quadrant signal for each quadrant is passed to a fan-in/fan-out unit (a different quarter of the unit for each quadrant) and for each quadrant, we use three copies of this signal.

- One goes directly to the trigger box where it is downscaled (downscaled particles).
- The second goes into a strobed coinc unit and from there to the trigger box. The strobed coinc unit gates each input (one for each quadrant of the CD) with the gamma gate (see section 4). This gives the particle- γ coincidence.
- The third copy goes to a scaler.

Scaler 69 has the delayed trigger for quadrant 1 and scaler 77 has the particle- γ coincidence for that quadrant. Scalers 70 and 78 have the corresponding signals for the second quadrant etc.




Figure 3: The generation of the triggers. The gamma gate is common to all four quadrants of the CD, but the rest is repeated for each quadrant with the second quadrant using channels 3 and 4 of the trigger box, the third using 5 and 6 and the last quadrant using the two remaining channels 7 and 8.

6 Generation of the DAQ dead signals




Figure 4: The generation DAQ dead signals.

This bit was changed quite a lot between July and November 2004. We use seven DAQ dead signals, three of which are inverted. These come from the top two halves of the fan-in/fan-out in R3.C3.S23. The four normal outputs are used to veto the LeCroy 465 AND modules in R3.C2.S13 and R3.C2.S15 (one for each quadrant of the CD - note that the ANDs are triple modules, so the three quadrants are on one module and the fourth on the other). The three inverted outputs are used to generate the “not busy” signal (R4.C4.S11), the “DAQ trig if DAQ not busy” (R3.C2.S13) and the “Ion and not DAQ dead” (R3.C4.S7).

7 Generation of the DAQ trigger




Figure 5: The generation DAQ trigger. The ionisation chamber part is shown dashed because it was not present for all experiments.

We use the DGF SYNCH to trigger readout (having ensured that we force a SYNCH at the end of the on/off window etc. - see section 19) However, if we send a DAQ trigger while the DAQ is busy, the trigger

module ignores it, so we need to postpone the DAQ trigger until the DAQ is no longer busy, otherwise it hangs. Since the SYNCH remains set until readout occurs, we AND this signal with the inverted DAQ dead signal, so that the moment we have both SYNCH and the DAQ is live, we generate a DAQ trigger.

The NIM \rightarrow ECL converter module is the same one used for the ADC busy lines and we have a special flat cable which connects the four ADC busy lines and the DAQ dead output from the VME trigger module to different channels of this unit.

It seems that between July and November 2004, this bit was changed a little, presumably to prevent deadlocks in the DAQ.

8 Generation of the ADC gates and the TDC stops




Figure 6: The generation of the ADC gates and TDC stops. Note that we have separate gates for each quadrant of the CD. This diagram shows the connections for the first quadrant. For the second quadrant we use channels 3 and 4 of the trigger box and so on.

In 2003 the ADC gate and TDC stop was the same for all parts of the CD, but in 2004 we changed it to have separate gates for each quadrant.

We want to generate an ADC gate, a TDC stop and a signal for the timestamper DGF corresponding to a given quadrant (four separate gates, four separate stops and four timestamping DGFs) all of the following conditions are satisfied:

- The trigger box gave a signal for that quadrant. Note that we can have both downscaled particle and particle- γ coincidence for the quadrant and these need to be ORed together.
- We are in the on window or the off window. This is the equivalent of the GFLT signal sent to the DGF.
- The ADC for that quadrant is **not** busy.
- The DGF is **not** busy (indicated by the DGF SYNCH signal).
- The DAQ is **not** dead (indicated by the VME trigger module).

The ADC busy signal is provided via a special flat cable which takes the ECL signal from each ADC to a single ECL \rightarrow NIM converter. Note that the DAQ dead signal coming from the VME trigger module is also on the same cable and uses the same ECL \rightarrow NIM converter.

Note also that in 2003 we generated an artificial ADC busy which started at the same time as the real one, but lasted longer. In 2004 we used the real ADC busy signal.

The ADC gate is taken directly from the fan-in/fan-out module at the end of figure 6 and sent to the CAEN V785 ADC gate input. Then another cable goes from the gate output to the TDC corresponding to the same quadrant.

9 Generation of the EBIS window

Note: The “EBIS window” is the window provided by a signal from the EBIS saying there is a bunch of particles incident on the target. The terms “on window” and “off window” refer to two measuring periods of the same length, one on beam the other off beam. Consequently “EBIS window” and “on window” are synonymous terms.




Figure 7: The generation of the EBIS window (on window)

A TTL signal is sent from the EBIS which after conversion to NIM is used to generate an 800 μs gate. The length of this gate needs to be set for the particular experiment to match the opening time of the EBIS gate.

10 Generation of the off window




Figure 8: The generation of the off window

We want to open one window when the EBIS gate is open (the on window) then read out the data and then open a second window (the off window) in the gap between EBIS pulses and read that data out in time to be ready for the next EBIS pulse. The on window and off window must be of equal length, so the off window length needs to be set to same as the EBIS gate.

When the EBIS signal comes, we send a 3 μs pulse to the timestamping DGF (note that this was 5 μs in 2003) and start a gate which is closed by the end of readout (the readout being forced at the end of the on window). We also start another gate which is closed before the next EBIS pulse (10 ms in this case, but again this has been 12.25 ms in another experiment). This is used to truncate the off window, so we should also allow some time for reading out the off window data, in order to be ready to acquire when the next EBIS pulse comes. However, normally, it should be possible for the off window to close and readout to complete long before the next EBIS pulse.

In the setup, we use the fact that the end of the on window triggers readout, causing the DAQ to become dead for a while, so we wait until the DAQ is not dead again after the on window.

We start the off window 60 μs after the end of readout (obtained from the disappearance of the DAQ dead signal) to allow the ADCs and TDCs time to start. Note that originally, 1 μs but it seems that this is not enough time for everything to be ready. Since we have 10 ms in which to open the window, this extra delay isn't a problem. Until we increased this value, the DAQ would hang from time to time. Afterwards it didn't hang nearly as much.

The on window and off window should be of equal length. It is possible to have an on window without an off window, but not the other way round, since the off window is triggered by the first readout within about 10 ms after the on window. Consequently, if there is no on window, the off window also disappears.




Figure 9: The generation of the on/off window

11 Generation of the on/off window

The on/off window is simply the logical OR of the on window and the off window. It is used to generate the GFLT, so the terms “on/off window” and “GFLT” are synonymous.

12 The generation of the GFLT




Figure 10: The generation of the GFLT from the on/off window

The global first level trigger (GFLT) signal which is sent to each DGF is the same thing as the on/off window. We use two 3 fan-in 39 fan-out modules to distribute the signal.

Typical times would be 800 μ s on window then about 3.2 ms readout and finally 800 μ s off window. The on and off window lengths should be equal and fixed. The readout time varies depending on the events.

13 The DGF BUSY/SYNCH loop




Figure 11: The DGF BUSY/SYNCH loop

The DGF BUSY/SYNCH loop is made using two 39-fan-in/3-fan-out modules and two 3-fan-in/39-fan-out modules. We take all the BUSY outputs from the DGFs and feed them into the fan-in inputs and then take the outputs to the fan-outs and send the result to each SYNCH input. In this way, if any one DGF is BUSY, the SYNCH line is set to logic one and if all the DGFs are acquiring it is logic zero.

14 Note about special signals

We get several signals from different parts of the REX ISOLDE setup. When ISOLDE is running, we get a T1 signal when a bunch of protons hits the ISOLDE target, then a T2 signal when the gate is opened allowing beam from ISOLDE into REXTRAP. Since these two signals are correlated one-to-one with a fixed time delay between them, it is not necessary to have both. If the laser is running in on/off mode, we get a signal indicating laser on or laser off (sent into the pattern unit).

Note: The T1 and T2 signals are different for each separator (HRS and GPS). You need to hook them up in the ISOLDE control room for the appropriate separator.

The individual bunches of protons hitting the ISOLDE target form part of a supercycle and we receive the PS signal indicating the start of each supercycle. Comparing the T1 pulses to the PS signal indicates how many bunches we receive in each supercycle.

Once the ions have been trapped in the REXTRAP and ionised in the EBIS, the EBIS gate is opened to allow them into the REX accelerator. We receive the EBIS gate signal which indicates when ions are released from the EBIS and which is directly correlated to our beam. Consequently, we use this signal to generate the on-beam window.

15 Generation of the T1 timestamp

We want two signals:

- the timestamp of the T1 signal itself, which we obtain by converting it from TTL to NIM and then generating a gate which we send via a linear fan-in/fan-out to the DGF but vetoed by the DGF busy.
- The timestamp of a time $10 \mu\text{s}$ after the DGF starts acquiring after a T1 pulse.

In order to distinguish the two signals, they are combined using a linear fan-in/fan-out, but with one of them attenuated using a LEMO T-piece with a 50Ω terminator.

The T1 timestamp is produced in a similar way to the proton supercycle timestamp. Both use the “ $10 \mu\text{s}$ after DGF not busy” signal, so the bottom three modules in figures 13 and 12 are the same.




Figure 12: The generation of the T1 timestamp

16 Generation of the proton supercycle timestamp




Figure 13: The generation of the proton supercycle timestamp

We want to generate a DGF timestamp giving the time of arrival of the proton supercycle start signal and a time $10 \mu\text{s}$ after the first readout after a proton supercycle start signal.

We receive an TTL signal each time that we get a proton supercycle start which we feed into an TTL \rightarrow NIM converter. We do two things with this signal.

- Firstly, we use the signal to start a direct gate to indicate the start of the proton supercycle, which we send via a linear fan-in/fan-out to a timestamping DGF.
- Secondly, we want the timestamp for the first time the DGF is not busy after a supercycle start. To get this, we start a timer when the supercycle starts and the DGF is busy, and reset that timer $10 \mu\text{s}$ after the DGF is no longer busy. Then we use the end marker of that timer to generate a gate, which is sent to the same timestamping DGF via the linear fan-in/fan-out.

Note that the bottom three modules figures 13 and 12 are the same.

Why don't we attenuate one of the two inputs to the linear fan-in with a 50Ω piece in order to get different pulse heights for the two kinds of signal? We do this with the T1 signal and it seems that we should do it to PS as well. I am told this isn't really that important for PS.

I notice that in May 2005, the inverted output from the TTL to NIM convertor is being used. This must have already been like that in November 2004, but I'm not sure if it was like that in July 2004. My notes from July 2004 indicate we used the normal output, but this could be an error. In any case, now we use the inverted output.




Figure 14: The generation of the control signals for the scalers

17 Generation of the control signals for the scalers

We want to read out the scalers every second. We convert a 1 Hz clock pulse into NIM and use it to start a gate which we fan out to the control inputs for the three scalers.

However, we only want to read when the EBIS pulse comes (on window), so we use the end marker of our gate and reset it with the EBIS pulse.

18 Generation of the control signals for the pattern unit




Figure 15: The generation of the control signals for the pattern unit. The four inputs come from the LeCroy 365 AND modules which generate the signal to produce the ADC gate for each quadrant, which is already vetoed by DAQ dead, ADC busy etc.

The pattern unit needs to record the pattern when each ADC gate occurs. However, we only have one pattern unit and four separate ADC gates (one for each quadrant) so we need to OR these signals logically. We take the signal out of the AND in the middle of figure 6 for each quadrant and pass them into a fan-in/fan-out. The resulting signal is the control for the pattern unit.

There was another change here between July and November 2004. An extra AND module was inserted between the fan-in/fan-out module called “EBIS pulse” and the reset for the “1 s in EBIS” dual timer. The other part of the and comes from the output of that dual timer, so the timer is only reset if it is actually running.

19 Generation of the forced readout




Figure 16: The generation of the forced readout

We take the GFLT signal, delay it 50 μ s and then use it to generate a 10 μ s wide pulse which we send to the fan-in of the DGF BUSY/SYNCH loop (see section 13). This causes the DGFs to finish their run and since the SYNCH is used to generate the DAQ trigger (see section 7).

20 The pattern unit bits

The first eight bits of the pattern unit correspond to the eight channels of the trigger box. i.e. downscaled quadrant 1, quadrant 1 and γ , downscaled quadrant 2 . . .

In August 2004, we added in the laser on/off signal (bit 9) into the pattern unit. Note that this bit indicates the request that we send to the laser cabin (which is ignored in manual mode) not the actual state of the shutter.

Note that the highest bits should not be used as they are used internally in the acquisition program.

Bit	Signal
1	downscaled quadrant 1
2	quadrant 1 particle γ
3	downscaled quadrant 2
4	quadrant 2 particle γ
5	downscaled quadrant 3
6	quadrant 3 particle γ
7	downscaled quadrant 4
8	quadrant 4 particle γ
9	Laser on

21 Laser on/off

The laser on/off signal is generated from the PS signal, because we want to synchronise with a proton supercycle. The idea is to switch alternately one supercycle laser on, one supercycle with laser off.

In 2003 and 2004, a TTL signal was sent to the laser cabin and used to switch the laser on and off. The box in the laser cabin had a switch “TTL/MAN” which in TTL mode responded to this signal and in MAN mode remained either continuously open or closed according to the state of a second switch (which we left always “open”). The same signal which was sent as TTL to the laser cabin was sent as ECL to the pattern unit. The problem with this is that this means the acquisition records the state that the laser shutter should be in if set to TTL mode. However, if the switch is set to manual, the pattern unit still sees the alternating on/off cycle, even though the shutter is not moving. Furthermore, in 2004 we had the situation where the shutter itself jammed and then the acquisition recorded the state the shutter should have been in, not the actual state.

Another shortcoming of the setup in 2004 was the use of a timer to determine the length of time the shutter should stay open or closed. In fact, the supercycle length varies depending on the number of pulses per supercycle. Also, we don’t know if the time between the arrival of the PS signal and the opening of the shutter is short enough if we get the first pulse in the supercycle. It is possible that the laser shutter isn’t open fast enough, so that the first pulse is effectively in laser off mode, even though the bit says we are in laser on mode. This is not a problem if somebody else gets the first pulse and we get later ones.

In May 2005 a new system was installed. Now the PS signal is sent directly to some electronics in the laser cabin which does the switching automatically. This electronics is controlled by a control box in R2.C3.S23, to which it is connected by three cables. Depending on the position of the control on this module, the electronics in the laser cabin either sets “laser on”, “laser off” or “laser on/off” (i.e. switching at each supercycle) modes. This means we don’t have to go to the laser cabin in order to switch between laser on and laser on/off modes.

There is one additional TTL signal which goes through the TTL to NIM convertor in R3.C4.S1 and from there up to the pattern unit.

22 1 MHz and EBIS, 1 MHz and GFLT

For some reason these signals weren’t documented in 2004. We simply take the logical AND of the 1 MHz clock signal (after conversion from TTL to NIM) and the EBIS signal. Since the EBIS signal is equivalent to the on window, this counts the number of microseconds in the on window. This result is sent to scaler 91. Similarly, we take the AND of the 1 MHz clock and the GFLT, which is equivalent to the on/off window and send it to scaler 92.




Figure 17: The 1 MHz signals.

This means that we use scaler 90 to measure the time between scaler readouts in microseconds, scaler 91 to give the amount of time the on window was active in that period and scaler 92 gives the sum of the on window and the off window.

23 Ionisation chamber




Figure 18: The electronics for the ionisation chamber.

In August 2004, the Munich people installed an ionisation chamber on another beamline sharing electronics with the rest of the Miniball DAQ. Its electronics is fairly simple. A signal is fed into a TFA, then a CFD and then ANDed with the inverted DAQ dead (i.e. DAQ not dead). This signal is passed into a gate generator and used to gate an ADC. The energy signals are sent into linear amplifiers and then into that ADC via a LEMO to flat cable converter.

We take the inverted out signal from the gate generator used to generate the ADC gate and send it to the coincidence unit in R3.C2.S13 labelled “DAQ trig if DAQ not busy” and and it with the other signals,

removing the condition on the DGF SYNCH. This change can be made by pressing the buttons on the module. We remove the condition on the SYNCH and add the one on the ionisation chamber when we want to use the ionisation chamber and the other way round for normal operation. This is shown by a dashed line in figure 5.

Note that this electronics is used with a separate version of the DAQ, but it is set so that we just have to switch from one program to another without changing the electronics, other than pressing the two buttons on the coincidence unit mentioned in the previous paragraph.

In May 2005, a minor change was made, moving the ionisation chamber electronics from R2.C3 to R1.C3.

24 Scalers

There are three 32-channel scaler modules in use. Scalers 1 to 64 are for the PPAC. The third has a variety of different signals from different sources.

Scaler	Signal	Scaler	Signal
65	Q1 free	73	Q1 accepted
66	Q2 free	74	Q2 accepted
67	Q3 free	75	Q3 accepted
68	Q4 free	76	Q4 accepted
69	Q1 delayed	77	Q1 and gamma
70	Q2 delayed	78	Q2 and gamma
71	Q3 delayed	79	Q3 and gamma
72	Q4 delayed	80	Q4 and gamma
Scaler	Signal	Scaler	Signal
81	Q1 gate	89	GFLT
82	Q2 gate	90	1 MHz
83	Q3 gate	91	1 MHz and on win
84	Q4 gate	92	1 MHz and GFLT
85	EBIS pulse	93	T1
86	Total DGF	94	PS
87	Si OR	95	
88	SYNCH	96	

25 Positions of modules in crates and racks

The racks are numbered from 1 to 5 with 1 being closest to the wall. The crates are numbered from 1 with crate 1 being highest. The notation R1.C2 means rack one crate two and R1.C2.S15 refers to the module in slot 15 of that crate. These numbers are shown on the circuit diagrams.

25.1 Rack 1

In rack 1 we have the PPAC HV supply and one camac crate without a crate controller.

25.1.1 CAMAC crate R1.C2

- Slot 15 - LeCroy 4532 MALU (Quadrant 1 OR)
- Slot 17 - LeCroy 4532 MALU (Quadrant 2 OR)
- Slot 19 - LeCroy 4532 MALU (Quadrant 3 OR)
- Slot 21 - LeCroy 4532 MALU (Quadrant 4 OR)

25.1.2 NIM crate R1.C3

- Slot 15-16 CF4000 quad CFD “unused”, “unused”, “ionisation chamber”
- Slot 17-18 Ortec 454 TFA “ionisation chamber”
- Slot 19-20 TC243 linear amplifier “IC”
- Slot 21-22 TC243 linear amplifier “SI”
- Slot 23-24 Lemo to flat cable converter for ADC for ionisation chamber.

25.2 Rack 2

In rack 2 we have two crates for the RAL 109 modules with their power supplies and a low-power NIM crate.

25.2.1 KM-6 crate R2.C1

- Slot 1 - RAL 109 s/n 70 (annular strips)
- Slot 2 - RAL 109 s/n 71 (annular strips)
- Slot 3 - RAL 109 s/n 72 (annular strips)
- Slot 4 - RAL 109 s/n 73 (annular strips)
- Slot 5 - RAL 109 s/n 74 (annular strips)
- Slot 6 - RAL 109 s/n 75 (annular strips)
- Slot 7 - RAL 109 s/n 76 (annular strips)
- Slot 8 - RAL 109 s/n 77 (annular strips)
- Slot 12 - RAL 109 s/n 90 (PAD detector)
- Slot 13 - RAL 109 s/n 91 (unused)

25.2.2 KM-6 crate R2.C2

- Slot 1 - RAL 109 s/n 78 (sector strips)
- Slot 2 - RAL 109 s/n 79 (sector strips)
- Slot 3 - RAL 109 s/n 80 (sector strips)
- Slot 4 - RAL 109 s/n 81 (sector strips)
- Slot 5 - RAL 109 s/n 82 (sector strips)
- Slot 6 - RAL 109 s/n 83 (sector strips)
- Slot 7 - RAL 109 s/n 84 (sector strips)
- Slot 8 - RAL 109 s/n 85 (sector strips)
- Slot 9 - RAL 109 s/n 86 (sector strips)
- Slot 10 - RAL 109 s/n 87 (sector strips)
- Slot 11 - RAL 109 s/n 88 (sector strips)
- Slot 12 - RAL 109 s/n 89 (sector strips)

25.2.3 NIM crate R2.C3

- Slot 1-2 ECL \rightarrow NIM converter
- Slot 3-4 Silena quad bias supply (CD bias)
- Slot 5-6 Silena quad bias supply (CD E detector bias (not used))
- Slot 7-8 Ortec 480 pulser (CD pulser)
- Slot 9-10 12 Volt power
- Slot 11-12 Triple preamp power supply “beam dump detector”, “ ΔE 1”, “ ΔE 2”
- Slot 13-14 Canberra dual HV supply (beam dump detector)
- Slot 15-16 CF4000 quad CFD “unused”, “unused”, “ionisation chamber” “unused”
- Slot 17-18 Ortec 454 TFA “ionisation chamber”
- Slot 23-24 Laser control box

25.3 Rack 3

In rack 3 we have a VME crate and three high-power NIM crates.

25.3.1 VME crate R3.C1

- Slot 1 - power PC
- Slot 2 - VME trigger module
- Slot 4 - CAEN V785 AG ADC
- Slot 5 - CAEN V785 AG ADC
- Slot 6 - CAEN V785 AG ADC
- Slot 7 - CAEN V785 AG ADC
- Slot 8 - CAEN V785 AG ADC (for ionisation chamber)
- Slot 9 - CAEN V775 TDC
- Slot 10 - CAEN V775 TDC
- Slot 11 - CAEN V775 TDC
- Slot 12 - CAEN V775 TDC
- Slot 13 - pattern unit
- Slot 14 - scaler with ECL inputs (PPAC Y)
- Slot 15 - scaler with NIM inputs
- Slot 16 - scaler with ECL inputs (PPAC X)
- Slot 17 - Wiener VC32
- Slot 18 - Wiener VC32
- Slot 19 - Wiener VC32
- Slot 20 - Wiener VC32 (not used)
- Slot 21 - VDIS

25.3.2 NIM crate R3.C2

- Slot 1-6 - TB8000 trigger box (triple width) “Trigger box”
- Slot 7-8 - CAEN N454 4x4 fan-in/fan-out “Delay trg Q1...4”
- Slot 9-10 - LeCroy 370C strobed coinc “coinc PG1...4”, “four channels unused”
- Slot 11-12 - CAEN N455 quad coinc logic “Q1...4 coinc or downscaled”
- Slot 13-14 - LeCroy 465 triple coinc unit “unused”, “DAQ trig if DAQ not busy”, “DGF Synch or ADC busy for Q4”
- Slot 15-16 - LeCroy 465 triple coinc unit “DGF Synch or ADC busy for Q1...3”
- Slot 17-18 - GG8000 octal gate generator “Q1...4 ADC/TDC gate”, “Ion chamber”, “Bad”, “DAQ trigger”, “ γ OR”
- Slot 19-20 - LeCroy 429A 4x4 fan-in/fan-out “Q1...4 ADC/TDC gate”
- Slot 21-22 - IKP modified TFA (in Ortec housing) “Mult out”
- Slot 23-24 - CAEN N454 2x8 fan-in/fan-out “Gamma Or”, “Si OR”

25.3.3 NIM crate R3.C3

- Slot 1-2 - TTL/NIM → ECL converter (set to NIM) “DAQ Trigger”
- Slot 3-4 - ECL → TTL/NIM converter (set to NIM) “DAQ dead”, “ADC dead”
- Slot 5-6 - LeCroy 429A 4x4 fan-in/fan-out “ADC 1 busy”, “ADC 2 busy”, “ADC 3 busy”, “ADC 4 busy”
- Slot 7-8 - Noname dual timer “Delay Q1 (both parts)”
- Slot 9-10 - CAEN N93B dual timer “Delay Q2 (both parts)”
- Slot 11-12 - Noname dual timer “Delay Q3 (both parts)”
- Slot 13-14 - CAEN N93B dual timer “Delay Q4 (both parts)”
- Slot 15-16 - LeCroy 429A 4x4 fan-in/fan-out “Q1...4 trigger”
- Slot 17-18 - ECL → NIM converter “TDC starts”
- Slot 19-20 - LeCroy 365 AL four-fold logic unit “unused”, “1s and EBIS”
- Slot 21-22 - CAEN 2255A Dual Timer “?”, “unused”
- Slot 23-24 - CAEN N454 4x4 fan-in/fan-out “DAQ dead”, “1s in EBIS”

25.3.4 NIM crate R3.C4

- Slot 1-2 - TTL → NIM and NIM → TTL converter “EBIS pulse”, “PS”, “T1”
- Slot 3-4 - Clock “1 Hz”, “1 MHz”
- Slot 5-6 - LeCroy 429A 4x4 fan-in/fan-out “not used”, “1 MHz”, “not used”, “not used”
- Slot 7-8 - CAEN N455 quad coinc logic “Ion and not DAQ dead”, “1 MHz and on win”, “DGF busy and PS start”, “1 MHz and GFLT”
- Slot 9-10 - Noname dual timer “Delayed TS PS start”, “PS del”
- Slot 11-12 - CAEN 2255B dual timer “unused”, “direct TS PS start”
- Slot 13-14 - CAEN N93B dual timer “1 sec in EBIS”, “TS gate T1 delay”
- Slot 15-16 - CAEN N454 4x4 fan-in/fan-out “T1”, “DGF busy”, “10 μ s after DGF not busy”, “start PS cycle”
- Slot 17-18 - CAEN N93B dual timer “T1 veto”, “TS delay”
- Slot 19-20 - CAEN 2255B dual timer “10 μ s del”, “DAQ trigger”
- Slot 21-22 - Philips quad linear fan-in/fan-out “PS TS”, “T1 not busy”
- Slot 23-24 - CAEN 2255B dual timer “TS gate T1”, “Force readout”

25.4 Rack 4

In rack 4 we have three CAMAC crates and one high-power NIM crate.

25.4.1 CAMAC crate R4.C1

- Slot 4 - XIA DGF 1159
- Slot 5 - XIA DGF 1171
- Slot 6 - XIA DGF 1139
- Slot 7 - XIA DGF 1161
- Slot 8 - XIA DGF 1162
- Slot 9 - XIA DGF 1163
- Slot 10 - IKP 40 MHz clock

- Slot 11- XIA DGF 1107
- Slot 12 - XIA DGF 1101
- Slot 13 - XIA DGF 1100
- Slot 14 - XIA DGF 1104
- Slot 15 - XIA DGF 1106
- Slot 16 - XIA DGF 1147
- Slot 17 - XIA DGF 1145
- Slot 18 - XIA DGF 1149
- Slot 19 - XIA DGF 1137
- Slot 20 - XIA DGF 1134
- Slot 21 - XIA DGF 1129
- Slot 22 - XIA DGF 1151
- Slot 24-25 Wiener CC32

25.4.2 CAMAC crate R4.C2

- Slot 4 - XIA DGF 1192
- Slot 5 - XIA DGF 1138
- Slot 6 - XIA DGF 1123
- Slot 7 - XIA DGF 1130
- Slot 8 - XIA DGF 1132
- Slot 9 - XIA DGF 1109
- Slot 10 - IKP 40 MHz clock
- Slot 11 - XIA DGF 1115
- Slot 12 - XIA DGF 1150
- Slot 13 - XIA DGF 1193
- Slot 14 - XIA DGF 1169
- Slot 15 - XIA DGF 1113
- Slot 16 - XIA DGF 1189
- Slot 17 - XIA DGF 1170
- Slot 18 - XIA DGF 1190
- Slot 19 - XIA DGF 1118
- Slot 20 - XIA DGF 1108
- Slot 21 - XIA DGF 1124
- Slot 22 - XIA DGF 1181
- Slot 24-25 Wiener CC32

25.4.3 R4.C3

- Slot 4 - XIA DGF 1154 “Q1 timestamp”
- Slot 5 - XIA DGF 1152 “Q2 timestamp”
- Slot 6 - XIA DGF 1142 “Q3 timestamp”
- Slot 7 - XIA DGF 1167 “Q4 timestamp”
- Slot 8 - XIA DGF 1176 “EBIS”, “T1”, “PS”
- Slot 9 - XIA DGF 1178 “Beam dump”
- Slot 10 - IKP 40 MHz clock
- Slot 11 - XIA DGF 1174
- Slot 12 - XIA DGF 1175
- Slot 13 - XIA DGF 1122
- Slot 14 - XIA DGF 1120
- Slot 15 - XIA DGF 1194
- Slot 16 - XIA DGF 1166
- Slot 17 - XIA DGF 1184
- Slot 18 - XIA DGF 1103
- Slot 19 - XIA DGF 1148
- Slot 20 - XIA DGF 1119
- Slot 21 - XIA DGF 1158
- Slot 22 - XIA DGF 1153
- Slot 24-25 Wiener CC32

25.4.4 NIM crate R4.C4

- 1-2 - Quad bias supply (ΔE detectors)
- 3-4 - CAEN N93B dual timer “EBIS”, “Max on/off window”
- 5-6 - CAEN N454 4x4 fan-in/fan-out “EBIS pulse”, “GFLT”, “not used”, “not used”
- 7-8 - LeCroy 365AL 4-fold logic unit “TDGF busy”, “End RDO and in time”
- 9-10 - CAEN N93B dual timer “TDGF Ebis”, “Off window”
- 11-12 - CAEN 2255B dual timer “Not busy”, “End RDO”
- 13-14 - IKP 3 fan-in/39 fan-out “GFLT”
- 15-16 - IKP 3 fan-in/39 fan-out “GFLT”
- 17-18 - IKP 3 fan-in/39 fan-out “Synch”
- 19-20 - IKP 3 fan-in/39 fan-out “Synch”
- 21-22 - IKP 39 fan-in/3 fan-out “Busy”
- 23-24 - IKP 39 fan-in/3 fan-out “Busy”

25.5 Rack 5

Rack five contains (from top to bottom) the filling computer monitor, the filling computer, its keyboard, the PT100 box, the four manifolds controller boxes (with box D highest), DAQ computer, raid array, UPS and CAEN HV power supply.