

Doc For Ionisation Chamber



At the beam line neighbouring to the MB beam line a dE-E-detector is mounted which allows to determine the composition of the beam. The detector consists of an ionisation chamber IC (dE) and a Si detector (E). The ionisation chamber is operated with P10 (90% Argon and 10% Methan). In front of the detector a attenuator reduces the beam intensity roughly by a factor 100.

Contents

Switch ON
DAQ
Operation
Switch OFF

The IC Si telescope is now mounted at the 0 degree beamline and equipped with a pressure regulation system. It is now operated with CF4 gas. The older documentation on cabling and DAQ given below remains valid.

Documentation of the pressure regulation system can be found here.

Basic operation instructions

The control unit (NIM module) which controls the piezo valve has two switches:

upper switch

- IST: the actual difference between the pressure in the IC and the reference pressure (prevacuum between prepump of turbo and turbo)
- SOLL: the target pressure value which the control unit should keep

lower switch

- ZU: piezo valve closed, no gas
- V: piezo valve open

!!!!!!! ATTENTION !!!!!!!!

The decimal point of the manometer display of the control unit is misleadingly shifted one position to the left, i.e. if you read 20.0 on the display this means 200 mbar!!!

!!

Evacuating

- lower switch closed (ZU) or gas bottle closed
- upper switch to IST to control pressure
- yellow valve between prepump and IC closed
- switch on prepump
- open slowly yellow valve

Operation

- prepump running, yellow valve open, lower switch is probably on ZU (closed)
- upper switch to SOLL
- choose the pressure you wish
- upper switch to IST
- yellow valve to value which corresponds to the pressure to be stabilised (see table)
- open gas bottle (1 bar should be sufficient)
- lower switch to V (open)


You'll now see that the control unit stabilises (after some over-/undershoot ...) the chosen pressure ...

End of operation (system remains ready to resume operation)

- high voltage off

- decrease slowly the HV of the Si down to zero
- decrease the HV of the IC down to zero
- switch both channels of HV unit to "Stand by"
- close the valve V1 at the bottle; valve V2 can remain open
- control that the pressure shown on P1 and P2 goes down
- you may now close the valve V6 of the beam line
- switch back the switches in the coincidence unit (A to ON, B to ON, and C to OFF)

After finishing return back to the normal directory of your current experiment, usually named cern-yymmdd, and start the normal DAQ.



Picture uses Nigel's convention of modules numbering: R==rack, C==crate, S==slot

Last edited on July 9, 2005 2:40 pm.